




Powering rare disease research and drug development with PacBio long-read sequencing

Nina Gonzaludo, PhD | Global Lead – Rare Disease
ngonzaludo@pacb.com

Who is PacBio?



Dr. Jonas Korlach



Dr. Stephen Turner


Nanofluidic design, single-molecule real-time chemistry
Cornell University, 1997




Menlo Park, CA




What is long-read sequencing?





SHORT READS
Missing sequence data
leads to gaps in genome
coverage and limits
variant detection



HiFi READS
Long reads map uniquely
and span large variants,
providing comprehensive
variant detection





What is HiFi sequencing?



HiFi read
>99.9% accuracy
15-20 kb (up to 100x longer than short reads)

Putting the “W” back in WGS



 Not all genome sequence data is created equally

Short-read WGS


-
-
-
-
- SNPs / small indels

- Needs a reference genome
- Miss 100s of millions of base pairs
- Blind to ~400 medically relevant genes in dark regions

HiFi WGS


- Structural variation
- Repeat expansions
- Segmental duplications
- Phasing/haplotype
- SNPs / small indels

- Clinical-grade, reference quality
- Resolve large, complex regions
- Methylation included

 PacBio HiFi delivers a new class of WGS

HiFi sequencing drives rare disease solves

What is the causal variant underlying disease? What is the disease mechanism of action?



Karyotype	Microarrays	Short reads		Long reads
		Exome	Genome	HiFi Genome
Chromosomal abnormalities	Copy-number variants >50kb	SNVs & indels, some large exonic variants	SNVs, indels, some large variants	SNVs, indels, SVs, CNVs, phasing, translocations, inversions, repeat expansions
~5% rare disease explanation rate	~10%	~30%	~40%	>50%

Familial long-read sequencing increases yield of de novo mutations

Michelle D Noyes¹, William T Harvey¹, David Porubsky¹, Arvis Sulovari¹, Ruiyang Li¹, Nicholas R Rose¹, Peter A Audano¹, Katherine M Munson¹, Alexandra P Lewis¹, Kendra Hoekzema¹, Tuomo Mantere², Tina A Graves-Lindsay³, Ashley D Sanders⁴, Sara Goodwin⁵, Melissa Kramer⁵, Younes Mokrab⁶, Michael C Zody⁷, Alexander Hoischen⁸, Jan O Korbel⁴, W Richard McCombie⁵, Evan E Eichler⁹

>25% more de novo SNPs and indels discovered with HiFi compared to short read WGS

ARTICLE · Volume 24, Issue 6, P1336-1348, June 2022 · [Open Access](#) [Download Full Issue](#)

Genomic answers for children: Dynamic analyses of >1000 pediatric rare disease genomes

Ana S.A. Cohen^{1,2,3}, Emily G. Farrow^{3,4}, Ahmed T. Abdelmouly⁴, Joseph T. Alaimo^{2,3}, Shivarajan M. Amudhavalli^{3,5}, John T. Anderson⁶, Lalit Bansal⁴, Lauren Bartik^{3,5}, Primo Baybayan⁷, Bradley Belden¹, Courtney D. Berrios¹, Rebecca L. Biswell¹, Pawel Buczkowicz⁸, Orion Buske⁸, Shreyasee Chakraborty⁷, Warren A. Cheung¹, Keith A. Coffman⁴, Ashley M. Cooper⁴, Laura A. Cross⁵, Tom Curran⁹, Thuy Tien T. Dang⁴, Mary M. Elfrink¹, Kendra L. Engleman³, Erin D. Fecke⁴, Cynthia Fieser⁴, Keeby Fitzgerald⁴, Emily A. Fleming⁵, Bendi N. Gaden⁵, Jennifer L. Cannon⁵, Rosa N. Gallego-Morales^{3,4}, Margaret Gibson¹

HiFi found 2x more structural variants, 4x more rare transmitted SVs compared to short read WGS

Other HiFi tools: Kinnex full-length isoform seq, methylation, chromatin accessibility

HiFi sequencing drives rare disease solves

NIH funds new *All of Us* Research Program genome center to test advanced sequencing tools

JUNE 14, 2023

PacBio and Radboud University Medical Center Team Up for New Discoveries in Rare Disease Research

FEBRUARY 29, 2024 | HUMAN GENETICS RESEARCH

Progress report: HiFi sequencing enables Care4Rare Consortium to advance rare disease research

Invitae and Pacific Biosciences Collaborate to Develop Whole Genome Sequencing-Based Assays for Pediatric Epilepsy Diagnostics

Kansas City, October 06, 2023

Children's Mercy Kansas City First to Use 5-Base HiFi Genomic Sequencing in the Clinical Setting

AUGUST 7, 2023

PacBio and GeneDx Launch Research Collaboration with the University of Washington to Study Long-Read Whole Genome Sequencing for Increased Diagnostic Yield in Neonatal Care

HudsonAlpha researchers use highly accurate long-read sequencing technology to help diagnose rare disease



Radboud Universiteit



HiFi sequencing for drug development: ASO design



Long read sequencing enables:

- Precise mapping of disease-causing mutations with base-pair level resolution
- Accurate identification of splice sites and regulatory elements
- Better understanding of local sequence context and target sites
- Improved detection of off-target effects, potentially improving safety profile


The diversity of *SNCA* transcripts in neurons, and its impact on antisense oligonucleotide therapeutics

 James R. Evans,  Emil K. Gustavsson,  Ivan Doykov,  David Murphy,  Gurvir S.Virdi,  Joanne Lachica,  Alexander Röntgen,  Mhd Hussein Murtada,  Chun Wei Pang,  Hannah Macpherson,  Anna I. Wernick,  Christina E. Toomey,  Dilan Athauda,  Minee L. Choi,  John Hardy,  Nicholas W. Wood,  Michele Vendruscolo,  Kevin Mills,  Wendy Heywood,  Mina Ryten,  Sonia Gandhi

doi: <https://doi.org/10.1101/2024.05.30.596437>



HiFi sequencing for drug development: gene editing





HiFi sequencing for drug development: gene editing



“SMRT sequencing can accurately measure the repeat expansion [...] to the C9orf72 locus in what has been a traditionally hard to sequence genomic region.”

“This is of particular value to sizing and phasing the repeat expansion and determining changes to the gene locus after gene editing”



Salomonsson SE, Maltos AM, Gill K, Arogundade OA, Brown KA, Sachdev A, et al. Validated assays for the quantification of C9orf72 human pathology. Sci Rep. 2024;14:828.

HiFi sequencing for drug development

Long-read sequencing enhances the discovery and development processes with high-quality data for confidence in results



PacBio HiFi sequencing supports:

- Target identification
- Biomarker discovery
- Library screening
- Directed evolution
- Cell line characterization
- Plasmid sequencing
- Animal model research
- Genotoxicity assessment

“PacBio long-read sequencing was and is fundamental to discovering Rose's drug, as well as other potential treatments for kiddos with diseases like hers. It provides a map for us early on that allows us to design the most impactful drug that we can at unprecedented speeds.”

- Casey McPherson, To Cure a Rose Foundation / Chrysalis Genetics



Sponsored by Ultragenyx

PacBio

Thank You



Sponsored by Ultragenyx